

Base Molds

Save Time, Energy and Money



Tissue-Tek® NanoMolds™ produce uniform paraffin blocks much quicker and release easier than traditional methods – all without the use of messy chemical mold release solutions. NanoMold stainless steel base molds are manufactured with an exclusive and patented Nano Technology coating that speeds up the embedding process helping laboratories save time, labor and increase efficiency. Even more time savings are achieved because residual paraffin is no longer left in the molds, drastically reducing the amount of cleaning necessary. Tissue-Tek® NanoMolds™ are available in five popular sizes to fit all of your embedding needs.

- Faster Block Release
- Eliminates Need For Using Mold Release Solution
- Less Cleaning
- Available In 5 Popular Mold Sizes

Specifications

(Dimensions sizes in mm's)

#4215 NanoMold, 9 x 9 x 4

#4216 NanoMold, 15 x 15 x 4

#4217 NanoMold, 20 x 20 x 5

#4218 NanoMold, 20 x 26 x 5

#4219 NanoMold, 20 x 33 x 5



Trusted Innovation™

Sakura Finetek USA, Inc. 1750 214th Street Torrance, CA 90501 Phone: (800) 725-8723

none: (800) 725-872 www.sakuraus.com